

PPLN Cavity Mixers:

Cavity configuration is an alternative way to enhance nonlinear frequency conversion. For your application convenience, we have developed cavity mixer platform which can be customized to a variety of mixer applications, such as external pump OPO (EP-OPO), Intra-cavity OPO (IC-OPO), intra-cavity SFG (IC-SFG), intra-cavity DFG (IC-DFG) etc. Example applications are for generating NIR signal wavelengths between 1.4-2 μm and MIR idler wavelengths between 2.3-4.5 μm .



Key features

- Cavity enhanced & high-efficiency, optimized for your specified input pump or desired output signal/idler
- Available output wavelength from UV/Visible to NIR/MIR
- Available mixing configuration from fundamental type to advanced type (such as IC-OPO, IC-SFG, IC-SHG, IC-DFG, EP-OPO etc.)
- Available for fiber or free space as input/output coupling interfaces
- Available with integrated thermistor/TEC for QPM temperature optimization & optional photodiode (PD) for power monitoring/automatic power control
- Convenient, compact and robust and available for a variety of application customization

Example: IC-OPO

an ultra-low power consumed optical parametric oscillator (OPO)


The intracavity design of the OPO utilizes the high circulating power in the cavity to achieve extremely low threshold pump power. 4 series with 3 wavelength cover ranges are provided upon request.

OPO-B and OPO-S are the OPO module with fixed wavelength specified by the customer in the range of either α , β or γ . In contrast, OPO-TB and OPO-TS are the wavelength tunable OPO that cover the entire range.


- B: broad bandwidth (few nm)
- S: single longitudinal mode
- T: tunable
- e.g. OPO-TS: tunable, single longitudinal mode

Parameters	unit	OPO-B	OPO-TB	OPO-S	OPO-TS
Signal/Idler wavelength	nm		α series: 1560-1880/2500-3300 β series: 1495-1640/3000-3700 γ series: 1440-1510/3600-4100		
Signal/Idler output power ^{*1}	mW		α series: 250/100 β series: 250/90 γ series: 200/70	^{*3} α series: 150/80 β series: 150/70 γ series: 100/40	
linewidth ^{*2}	GHz	300	300	<0.2	<0.2
Frequency stability	MHz/hr	N/A	N/A	<400	<400
Beam Quality	-		TEM ₀₀ , signal M ² <1.2, idler M ² <1.5		
Beam divergence	mrاد		Collimated, <1 mrad		
Power stability (rms)	%		<5		
Polarization			Linear, >20dB		

^{*1}. Peak value, power may vary with wavelength
^{*2}. The linewidth changes a little with wavelength
^{*3}. Preliminary specifications



Frequency stability of the OPO-S- β



power scaling of OPO-B- α

power scaling of OPO-B- α

Example available configurations:

