Leading Manufacturer of Crystals and Optics for Laser Applications

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Thank you for using CASTECH's products.

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Fujian CASTECH Crystals INC. is located in Fuzhou city, Fujian Province, in China. It is devoted to the development, manufacturing and marketing of nonlinear optical crystals, laser crystals and laser optics. CASTECH's products are widely used in the fields of laser and optical communication. In Mar.19, 2008 CASTECH went public in The Shenzhen Stock Exchange (Stock Code 002222).

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After twenty years untiring effort, CASTECH has developed flux, Czochralski and water solution crystal growing technology. We have IAD, IBS and EB coating techniques to meet different requirements. A whole set of crystal inspection and manufacturing controlling equipments including Zygo, Nikon, Trioptics ensures that every piece of CASTECH's product are carefully inspected. In 2001, CASTECH got the ISO 9001 certification, and in July 2006, the EU RoHS (Restriction of Hazardous Substance Directive) standard (2005/95 / EC) was introduced to the manufacturing process, which ensure that all of the products meet the RoHS requirements.

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Over 80% of CASTECH's products are exported to the USA, Japan, Germany and other countries and districts, and are praised as "China Crystal". Today CASTECH will work towards the goal of becoming the best global supplier for non-linear optical crystal and laser crystals.

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Thank you again for your cooperation. We hope you will continue to support us. Please read all of the product instructions in detail before using our products.

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
Disclaimer Statement:

- (1) This User's Manual was prepared by CASTECH's Technicians according to the product's characteristic. We intend to help the users to use our products more convenient.
- (2) The users should have knowledge of laser and optics, and are familiar with the laser manufacturing and related product manufacturing.
- (3) This User's Manual doesn't include all of the notes that the users should pay attention to. The user should avoid the influence and hurt by electricity, water, noise, strong light, etc. during assemble process.
- (4) The users may need more strict control when you manufacture different product.
- (5) It doesn't mean unqualified product if the product you receive has different mark with this User's Manual.

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
Example of Crystal Product:


①-Batch No.	⑤-Coating
②-Products and quantity	⑥-ID
③-Orientation	⑦-Serial Number
④-Dimensions	⑧-Manufacturing Date

Meaning: LBO, cut angle is ($\Theta = \Phi =$), with coating AR-1064/532nm, Serial NO. is 1-12345-0325.

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
Example of Optical products:




①-Batch No.	⑤ -Coating
②-Products	⑥ -Mounting
③-Quantity	⑦- Manufacturing Date
④ -Dimensions	

Meaning: 1/2 waveplate with coating AR-526.5nm, Phi25.4mm mounting.

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


If you have already specified mark requirements in your orders, you may skip this section. If you haven't specified please refer to our regular standard as follows.

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Type	Description
<p>$T^* \geq 2.5\text{mm}$</p>	<p>“●→”: “●” + “→”</p> <p>“●”: The marked surface is perpendicular to the tuning axis.</p> <p>“→”Orientation: S1 to S2 (along the beam path) .</p> 
<p>$T < 2.5\text{mm}$</p>	<p>The blackened surface is perpendicular to the tuning axis.</p> 
<p>Serial number</p>	<p>For the thick crystal ($T \geq 4\text{mm}$), the serial number (eg.0025) will be marked on the side face as shown in the diagram.</p> 

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







*Surface X,Y,Z is perpendicular to the crystal's optical axis.

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<p>BBO with fused silica substrate</p>	<p>S1 :BBO S2: fused silica</p> <p>substrate</p> <p>BBO</p> <p>holder</p>
<p>Waveplate mounting</p>	<p>The platform is parallel to the fast axis</p> <p>platform parallel</p> <p>axis</p>
<p>Oven and holder</p>	<p>The marked surface faces the positioning and power lines.</p> <p>crystal</p> <p>dispensing hole</p> <p>holder</p> <p>Oven</p> <p>Positioning</p> <p>Power line</p>
<p>Mounting of glued Crystals (Nd:YVO4+KTP)</p>	<p>Nd:YVO₄ a surface faces the notch of the holder.</p> <p>Si</p> <p>KTP</p> <p>Nd:YVO₄</p>

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
Products Instruction ■ ■ ■


Once the inner plastic box has been opened, CASTECH's engineer recommends to provide a dry and clean environment.

Temperature: $23 \pm 2^{\circ}\text{C}$

Humidity: $\leq 55\%$

The above conditions are required for product assemble process. The conditions in your laser system should be more strict than those to ensure long lifetime of your laser system. Especially when using the deliquescence crystals such as LBO, BBO, DKDP etc, please ensure that the cavity environment is dry, sealed, clean and the humidity had better than 15%.

- ◆ It is recommended that only qualified personnel can open the inner plastic boxes. Finger prints, oil and other substances should not contact the polished or coated surfaces directly.
- ◆ Package opening instructions:
Please cut the plastic bags and remove the plastic packaging. Once this has been done, carefully open the box along the four corners of the plastic packaging, paying special attention to opening the lid to avoid the crystal (optics) falling (it may be adhere to the lid during transportation); at the same time do not force open the lid to avoid breaking the crystal (optics) inside.
- ◆ Instructions of clamping the crystals (optics):
Using carbon tweezers or wearing a fingerstall to clamp the crystal (optics) without touching the polished or coating surfaces. The use of metal tweezers is not recommended to avoid breaking or scratching the crystals (optics).
- ◆ Instructions of cleaning the crystals (optics):
If the surfaces are contaminated, an air ball is recommended to blow on the surfaces. If there is still contamination on the crystal (optics) surfaces, please clean the surfaces with cleaning liquid and soft silk. We recommend using a mixing liquid of 50% high purity alcohol and 50% high purity ether. If you can not clean the crystal(optics), please do not hesitate to contact us.

- ◆ We suggest that do not expose the crystals (optics) to the air for more than thirty minutes if possible. Please refer to STORAGE sections for more details.

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USE of Nonlinear Optical Crystals:

- 1) Please make sure your laser power density is under the damage threshold of crystals.
- 2) Please install the crystals according to the crystal phase-matching type and polarization directions of fundamental lasers, and carefully adjust the crystal along the tuning angle.
- 3) Please avoid directly touching the crystal surface. If the surface need cleaning, it should be done by professional staff according to the cleaning manual.
- 4) For crystals used in UV band, it is important not to use package containing silicon (Si) element.
- 5) A heating and cooling rate of less than 3 °C / min is needed to avoid the crystal or the coating film cracking.
- 6) Handling carefully during assembly, too strong pressure may cause the cracking of crystal.
- 7) Align the laser beam to the middle of the crystal. If the laser beam hits the side of the crystal, it may volatile something and cause pollution of the laser cavity and damage the crystal surface.

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Laser Crystal

1. Ensure that the pumping wavelength is consistent with the absorption wavelength of the laser crystal.
2. When installation birefringent laser crystal such as Nd: YVO₄ or Nd: GdVO₄, it is important to make sure the polarization direction of the pump source is along the maximum absorption direction of the laser crystal.

(For information about crystal heating and cooling rates, the handling of crystal and alignment of laser beam requirements, please refer to “nonlinear crystal” 5-7)

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Green and blue Laser Kit

1. Make sure that the wavelength of pumping light is consistent with the absorption wavelength of the laser crystal.
2. It is recommended that the pumping power should be under 500 mW for green laser kit.
3. Ensure that when the crystals are in use, they are well cooled.

(For information about crystal heating and cooling rates, the handling of crystal and alignment of laser beam requirements, please refer to “nonlinear crystal” 5-7)

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Storage

1. Crystals (optics) should be well cleaned after use and placed into the box. At the same time non-contact packaging is required to prevent scratches or smudges on the two light passing surfaces. Ensure that the packaging is sealed, neat, moisture-proof, shockproof and anti-static.
2. Crystals (optics) should be stored under dry condition. They can be kept in vacuum packaging with desiccant to prevent from humidity and affects its normal function.
3. The recommended storage conditions; temperature $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$, humidity $< 20\%$, clean (maintained at 1000), and isolated.

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Service

If you have any questions with testing and debugging, please contact us without hesitate and provide information including PO No., the batch no. and the serial number on the packaging box or on the side face. We will seek to provide professional solutions for you as soon as we can.

When you order special crystals or optics, please inform us with the main parameters of the laser, such as pulse energy, pulse width, repetition rate, continuous power, beam diameter, beam mode, beam divergence and tunable wavelength range.

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
Appendix

The introduction of standardized elements, spare parts, technical documents etc. is covered in the appendix.


WWW.CASTECH.COM
MM'CV2TECH'COM

LBO

BBO


KTP

Nd:YVO4

OPTICS

THIN FILM FILTERS

ADD: 155 Yangqiao West road, Fujian 350002 P.R. China.

TEL: +86-591-83710533

FAX: +86-591-83711593

E-Mail: sales@castech.com

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Fujian CASTECH Crystals, Inc.